## **Greenway Daily Equipment Check List**

	Daily Checks w/c	Monday	Tuesday	Wednesday	Thursday	Friday	Saturday	Sunday
1	<b>Drug Checks</b> Two people (1 must be a clinician), to check these drugs every day.	Daily	Daily	Daily	Daily	Daily	Daily	Daily
	Midazolam, Diazepam & Morphine (Oramorph), Codeine							
	Oxycodone, Morphine Sulphate							
2	Consulting Room Boxes	Daily	Daily	Daily	Daily	Daily	Daily	Daily
	Check contents against list on box, report to SM any missing items.							
3	Panic Alarms	Daily	Daily	Daily	Daily	Daily	Daily	Daily
	All in a small box in host cupboard, ensure they are back in the box end of shift							
4	Patients and downloads folder (on desktop)			Daily	Daily	Daily	Daily	Daily
	Please make sure its clear start and end of shift							
5	Check respiratory hood and HCIS cleaning box are present and complete	Daily	Daily	Daily	Daily	Daily	Daily	Daily
	Check contents against list on both boxes							
	Clean contents and box itself with Clinell wipe							
6	Resus bag - Location: Store cupboard in corridor	Daily	Daily	Daily	Daily	Daily	Daily	Daily
	Check contents as per contents sheet, replace any stock if needed. Check any OO	D						
	and log any imminent. Ensure green tick is showing which indicates Defib is							
7	working. Clean defib with Clinell wipe	Della	Dette	Dette	Dethy	Della	Detha	Dethy
/	Blood Taking Box	Daily	Daily	Daily	Daily	Daily	Daily	Daily
	Ensure missing items are replaced from stock, check expiry dates and replace as							
_	required	D 1	D 11				D 1	D 11
8	Sharps boxes	Daily	Daily	Daily	Daily	Daily	Daily	Daily
	Should not exceed the 3-month date, should be closed but not locked, any							
	overfilled/not closed boxes to be reported							
9	Manual Scripts	Daily	Daily	Daily	Daily	Daily	Daily	Daily
	Count and record manual prescription packs in Audit book							
	Count and record manual prescription packs in Audit book							
10	Emergency drugs - Location: Inside drugs cupboard in store room	Daily	Daily	Daily	Daily	Daily	Daily	Daily
		,	,	,		,	,	
	Ensure drugs are in the appropriate drugs cupboard in store room							
11	Oxygen - Location: Small and Large bottles	Daily	Daily	Daily	Daily	Daily	Daily	Daily
	Please accurately indicate in the signature LARGE Cylinder Room 4							
	box the levels as below:							
	EMPTY / ¼ FULL / ½ FULL / ¾ FULL / FULL							
	Large Cylinder Driver's cupboard	u						
	Small Cylinder, Driver's cupboar	-d						
	Sindir Cymraer, Driver's Cupboar	~						

## **Greenway Daily Equipment Check List**

12	Medication issued from stock form									
	Ensure previous dated meds form is put in post box regardless of any meds issued									
13	Monitoring Drug Temperatures – Location: Drug Cabinet	Daily	Daily	Daily	Daily	Daily	Da	aily	Daily	
	Check temp and record Min and Max temp	MIN:	MIN:	MIN:	MIN:	MIN:	MIN:		MIN:	
	• If the alarms go off whilst on shift, please let the shift manager know	MAX:	MAX:	MAX:	MAX:	MAX:	MAX:		MAX:	
14	Handwash Audits	Daily	Daily	Daily	Daily	Daily	AM	PM	AM	PM
	Check which clinicians need handwash by logging into Radar Handwash Initial when completed handwash.									
	Weekly Checks	Monday	Tuesday	Wednesday	Thursday	Friday	Satu	ırday	Sunday	
15	Health and Safety Checks Weekly		Weekly							
	Complete H&S checklist, keep original in folder and flag any concerns									
16	Paediatric/Adult pulse oximeter - Location: Store Room within Sats monitor box Weekly	Weekly								
	Switch on machine, place on finger and ensure it takes a reading. Ensure Adult, child and Infant leads are available Clean the machine and box itself with Clinell wipe									
17	Nebuliser box			Weekly						
	Collect nebuliser box from the store room and ensure it is available for clinicians Check the machine powers on, and clean it with a Clinell wipe Ensure there is 1 x adult and 1 x paediatric nebuliser mask within the nebuliser box, checking the expiry dates and restocking where necessary									
18	Doppler Weekly					Weekly				
	Check Doppler is stored in cupboard and the equipment turns on. Clean the machine with Clinell wipe									
19	Blood Glucose Machine / Dual Keytone Weekly Check				weekly					
	Host to check machine as per guidelines Clean the machine and box itself with Clinell wipe									
20	Respiratory Hood needs charging once a week.				Weekly					
	Host to plug in the battery into charging unit for 3 hours See emailed instructions									
	BASE SPECIFIC CHECKS:	Monday	Tuesday	Wednesday	Thursday	Friday	Satu	ırday	Sun	day
21	Clinical Waste								We	ekly
	Please ensure the clinical waste bag is emptied if full, tagged and put in the clinical waste bin near staff entrance – keys on hook behind reception desk. This is essential before closing on a Sunday.									
	Sample safe checks	Daily	Daily	Daily	Daily	Daily	Da	aily	Da	nily
	Check of sample safe at the beginning and end of shift									

22	Headsets	Daily						
	Record how many and report any issues							
23	Panic Alarms- in reception Area and consulting rooms	Daily						
	Check Panic Alarm in reception area states ' System Ready'							

Week end date. .....